


INDIAN INSTITUTE OF TECHNOLOGY, BOMBAY

COURSE CH-117L	ROLL NO. 23B0912	NAME Aditya Sanapala	
ASSIGNMENT NO. 3		DUE DATE	SUB. DATE

Experiment -3: Determination of Phosphoric acid in Soft Drinks.

Aim: To calculate the amount of phosphoric acid in a soft-drink (Thumbrs up) using molybdenum blue.

Theory:



→ In this experiment, ammonium molybdate $(\text{NH}_4)_6\text{Mo}_7\text{O}_{24} \cdot 4\text{H}_2\text{O}$ forms a colorless hexavalent phosphomolybdate with the phosphate PO_4^{3-} ion.

→ In acidic conditions, the colorless complex is reduced to a blue pentavalent form by ascorbic acid.

→ The concentration of the pentavalent form can be determined spectrophotometrically by measuring its absorption at 830 nm.

Procedure:

1. Preparation of standard stock solution: 1.9412 g of KH_2PO_4 is dissolved in de-ionized water in a 1L volumetric flask and is topped off to 1L. (Standard A)

2. Preparation of working standard solution: Dilute the standard A 20 times by adding 2.5 mL A to 50 mL volumetric flask and top it to 50 mL (Standard B)

Dilute standard B 12.5 times by adding 4mL of standard B to

50mL volumetric flask and top it to 50mL (Standard C)

The concentration of standard C is 4ppm P_2O_5 .

3. Preparation of reducing solution R:

- Dissolve 1.00g of ammonium molybdate in 50mL de-ionized water (solution n)
- Dissolve 1.76g of ascorbic acid in 100mL de-ionized water (solution A)
- Dilute 17mL of $CC\cdot H_2SO_4$ in 200mL de-ionized water (solution S)
- Mix 39mL of solution M, 60mL of solution A, 125mL of solution S and top it upto 250mL in a volumetric flask.

4. Sample preparation (Solution X):

- Transfer 2.5mL of de-carbonated drink to a 50mL volumetric flask and top it up to 50mL.
- Prepare the sample test tubes as follows:

Vol (in mL)	Tube 1	Tube 2	Tube 3	Tube 4	Tube 5	Sample
Vol of Standard C	0	0.5	1	1.5	2	2.5
Vol of Soln X						0.25
Vol of Soln R	2	2	2	2	2	2
Vol. of water	3	2.5	2	1.5	1	0.5
Final Volume	5	5	5	5	5	5
Final C of P_2O_5 (in ppm)	0	0.4	0.8	1.2	1.6	2

- Heat these test tubes in a water bath for 10 mins.

INDIAN INSTITUTE OF TECHNOLOGY, BOMBAY

COURSE CH117L	ROLL NO. 23B0912	NAME Aditya Samapata	
ASSIGNMENT NO. 3		DUE DATE	SUB. DATE

→ Read the absorbance against a blank off at 830nm.

Observations:

- When N, A & S were added in order, there was a color change from colorless to yellow to pale green.
- After the test tubes were removed from the water bath, some of them had turned blue.
- Absorbance values for sample and standard tubes 1-6:

Solution	Concentration (ppm)	A ₈₃₀
Tube 1	0	0.0149
Tube 2	0.4	0.1148
Tube 3	0.8	0.1417
Tube 4	1.2	0.2065
Tube 5	1.6	0.2714
Tube 6	2	0.3838
Sample		0.2495

Calculations

- From the graph, equation of line is $y = 0.162n + 0.0149$
 - But, by simple linear regression, line is $\underline{y = (0.143n) n + 0.0693}$
- ∴ Concentration of sample is

$$0.2495 = (0.143n)n + 0.0693$$

$$\Rightarrow n = \frac{0.2495 - 0.0693}{0.1434}$$

$$\text{C. of } \text{B}_2\text{O}_5 = \boxed{n = 1.256 \text{ ppm}}$$

- Calculation of phosphorus content in sample:

Considering the dilutions made to the original sample and the transformation from B_2O_5 to P. The amount of P in the

$$\text{sample (in ppm)} = n \times 20 \times 20 \times \frac{62}{142}$$

$$= 1.256 \times 20 \times 20 \times \frac{62}{142}$$

$$\boxed{C = 219.357 \text{ ppm}}$$

Graph: Concentration vs. Age

23B0912 Aditya Samayala

$$y = (0.162)x + 0.0119$$

Liquid

$$x \square = 0.270$$

$$y \square = 0.35$$

