

COURSE CH117L	ROLL NO. 23B0912	NAME Aditya Samapala	
ASSIGNMENT NO. 2		DUE DATE	SUB. DATE <u>10</u> <u>10</u>

## Experiment-2: Thin Layer Chromatography

### Objectives:

- ① To prepare a fluorescent dye in microscale using a one-pot sequential amide formation - Nucleophilic aromatic substitution reaction.
- ② To understand the principles and application of Thin Layer Chromatography.
- ③ Analysis of TLC using 'ImageJ'
- ④ Determination of molecular weight by using TLC-MS.

### Introduction & Principles:

- Chromatography is a technique used to separate components of a mixture of compounds.
- Thin layer chromatography (TLC) is based on the selective adsorption of the solvent and solute on the stationary phase.
- The technique consists of applying the mixture on a stationary phase (a thin layer of silica gel,  $\text{SiO}_2$ , on aluminium) and allowing a mobile phase (toluene) to climb up the adsorbent.

- The most important criterium in qualitative identification is the  $R_f$  value (called retention value), which is defined as


$$R_f = \frac{\text{Distance traveled by the solute}}{\text{Distance traveled by the solvent}}$$

The dye:

- Fluorophores are molecules which can be excited and relax efficiently in a radiative pathway.
- On absorption of a photon, a molecule in the singlet ground state ( $S_0$ ) gets excited to the first excited state ( $S_1$ ).
- It can relax back to the ground state ( $S_0$ ) by emission of a photon of lower energy than the incident one.

# INDIAN INSTITUTE OF TECHNOLOGY, BOMBAY


COURSE CH171	ROLL NO. 23BPA12	NAME Aditya Sanapala	
ASSIGNMENT NO. 2		DUE DATE	SUB. DATE



Ground State

## Observation:

(under long IR)



- 4 spots are visible on the TLC when observed in long UV.
- From bottom, the spots are in:
  - light yellow
  - fluorescent (bright) yellow → This is our product.
  - greenish-yellow
  - and blue.
- The second spot from the bottom is our final product.
- After running the TLC through mass spectrometry, we obtain the molecular weight of the final product as 381.5 g/mol.
- After running the image of the TLC through "ImageJ" software, we get the area of the spots as follows:

Spot Color	Area
light yellow	21005.915
bright yellow (Product)	26814.643
greenish-yellow	10985.569
blue	15771.016

### Calculation:

Composition % of a particular spot =  $\frac{\text{integrated area of spot}}{\sum \text{integrated area of all spots}} \times 100$

# INDIAN INSTITUTE OF TECHNOLOGY, BOMBAY

COURSE CH1172	ROLL NO. 23B0912	NAME Abitya Samapala	
ASSIGNMENT NO. 2		DUE DATE	SUB. DATE

So, Composition % of pale yellow spot =  $\frac{21005.915}{74517.143} \times 100 = 28.18\%$

Composition % of bright yellow spot =  $\frac{26814.643}{74517.143} \times 100 = 35.98\%$

Composition % of greenish-yellow spot =  $\frac{10985.869}{74517.143} \times 100 = 14.74\%$

Composition % of Blue spot =  $\frac{15711.016}{74517.143} \times 100 = 21.08\%$

→ In mass spectrometry, radiation is NOT used but in NMR & IR Spectroscopy, radiation is used (NMR → radio, IR → IR).